



Port Caps

Create closed solvent delivery systems with customized multi-port closures. Available for many waste containers and tubing sizes.



7-port GL45 cap



3-port GL38 cap

Drum Feed Safety Clamps

A safe, fast, effective way to increase the safety of your chemical feed system. Fits standard plastic drums between 15 and 50 gallons



CFS-1117



Bench Top Solid Waste Container

Safely dispose of contaminated solid wastes (wet silica gel, kimwipes, pipette tips, gloves, and more) in a small waste container.



SWC-002R

SWC-002Y

TELLI

Audibly alerts to prevent spills and overflows. Fits most popular liquid waste containers.



CP-Telli-01



Talking Electronic
Liquid Level Indicator



eco funnel

WASTE DISPOSAL SOLUTIONS

for a safe & healthy laboratory

CPLabSafety



Features & Benefits

- Solves the “open waste container” problem
- Fits easily in a fume hood or sink
- Prevents fires by trapping volatile fumes
- Meets safety regulations
- Protects workers by minimizing exposure to odorless inhalants
- Widely used in academic, pharmaceutical, biotech, and government laboratories
- Considered a highly effective tool for chemical waste management by EH&S officers and industry professionals
- Choose a 4" ECO Funnel (450mL) or 8" ECO Funnel (2,000mL)
- 4 lid colors available to help prevent cross-contamination
- All models are customizable with HPLC waste line adapters
- Add up to 8 quick disconnect, luer lock, luer barb, or thread barb HPLC waste line adapters in various sizes



Think SAFETY



Gasket inside lid prevents 99.9% of chemical emissions

Large particle filter catches stir bars and other solids.

Small tube inside funnel releases pressure.

8" funnels come in 4 different lid colors



Various screw caps available to fit many common waste containers

Secondary tube on 8" models acts as a liquid-level indicator to prevent over-filling

Introducing the EF2, Our Newest Model

- 8" flat lid design deters accidental pouring while closed
- Simple latch design clicks open and shut easily with one hand
- Available for containers with 70mm and 83mm closures

