

Features of ECO Funnel®

The Simple, Safe and Effective Solution to Liquid Waste Management

Designed to make waste disposal simple and safe, ECO Funnel® safety features help prevent fires and spills in the lab



1. Gasket in the lid prevents 99.9% of emissions from escaping, reducing the risk of fire.
2. Variety of screw caps to fit common waste containers.
3. Removable filter separates large particles and stir bars.

4. The lid with latch makes closing the system convenient.
5. A small tube creates an air gap in the waste container to prevent over-filling.

ECO Funnels® Can Help You Comply with EPA Regulations:

29 CFR 1910.1450(e)(3)(ii): ...implement control measures to reduce employee exposure to hazardous chemicals

40 CFR 264.173(a): A container holding hazardous waste must always be closed during storage, except when it is necessary to add or remove waste

40 CFR 264.175(a)(b): Containers containing free liquids; A base must underlie the containers... The base must be sloped.... the containment system must have sufficient capacity to contain 110% of the volume of the containers or the volume of the largest container, whichever is greater...

